






Vindum MV Connector Model Listing

Model Number	Schematic	Connection Pattern	Dimensions (WxDxH) inch	MAWP (psi)	Tubing OD
MV-Union-HC		Straight Coupling	1.4"x1.4"x0.6"	12,000	1/8" or 1/16"
MV-TEE-HC		3-Way Coupling	1.4"x1.4"x0.6"	12,000	1/8" or 1/16"
MV-Cross-HC		4-Way Coupling	1.4"x1.4"x0.6"	12,000	1/8" or 1/16"
MV-5STA-HC (special order)		5-Way Manifold	1.4"x1.4"x0.85"	12,000	1/8" or 1/16"
MV-9STA-HC (special order)		9-Way Manifold	1.4"x1.4"x1.45"	12,000	1/8" or 1/16"